

Importing 3D Images

You've seen how to import a series of images as a 3D volume using the Volume Viewer app. You've also performed custom segmentations using the Volume Segmenter app. The Volume Viewer and Segmenter apps are useful tools to gain familiarity with your 3D images and perform common operations.

However, if you need to process many 3D images or perform highly customized operations, it might be more useful to analyze your volumes with code.

In this reading, you'll learn how to import, create, and view 3D objects in your scripts.

Table of Contents

Common Medical Image Formats.....	1
Reading Medical Images.....	2
3D Images as a Series of Images.....	4
3D Images as Volumes.....	5
Creating a Volume Object.....	6
Removing extra dimensions.....	6
Viewing a Volume Object.....	7

Common Medical Image Formats

Many medical scanners use file types specific to medical imaging. In addition to common image file formats, MATLAB also supports many medical imaging formats such as:

- DICOM ("Digital Imaging and Communications in Medicine")
- IMA ("Image", a disk image very similar to DICOM)
- NIfTI ("Neuroimaging Informatics Technology Initiative")
- Analyze 7.5

Let's take a closer look at the images we used in the "Analyzing 3D Images" video:

```
ls("./T2 MRI Scan")
```

```
.          T2_12.IMA  T2_142.IMA  T2_165.IMA  T2_188.IMA  T2_37.IMA  T2_6.IMA   T2_82.IMA
..         T2_120.IMA T2_143.IMA  T2_166.IMA  T2_189.IMA  T2_38.IMA  T2_60.IMA  T2_83.IMA
T2_1.IMA   T2_121.IMA  T2_144.IMA  T2_167.IMA  T2_19.IMA  T2_39.IMA  T2_61.IMA  T2_84.IMA
T2_10.IMA  T2_122.IMA  T2_145.IMA  T2_168.IMA  T2_190.IMA T2_4.IMA   T2_62.IMA  T2_85.IMA
T2_100.IMA T2_123.IMA  T2_146.IMA  T2_169.IMA  T2_191.IMA T2_40.IMA  T2_63.IMA  T2_86.IMA
T2_101.IMA T2_124.IMA  T2_147.IMA  T2_17.IMA   T2_192.IMA T2_41.IMA  T2_64.IMA  T2_87.IMA
T2_102.IMA T2_125.IMA  T2_148.IMA  T2_170.IMA  T2_2.IMA   T2_42.IMA  T2_65.IMA  T2_88.IMA
T2_103.IMA T2_126.IMA  T2_149.IMA  T2_171.IMA  T2_20.IMA  T2_43.IMA  T2_66.IMA  T2_89.IMA
T2_104.IMA T2_127.IMA  T2_15.IMA   T2_172.IMA  T2_21.IMA  T2_44.IMA  T2_67.IMA  T2_9.IMA
T2_105.IMA T2_128.IMA  T2_150.IMA  T2_173.IMA  T2_22.IMA  T2_45.IMA  T2_68.IMA  T2_90.IMA
T2_106.IMA T2_129.IMA  T2_151.IMA  T2_174.IMA  T2_23.IMA  T2_46.IMA  T2_69.IMA  T2_91.IMA
T2_107.IMA T2_13.IMA   T2_152.IMA  T2_175.IMA  T2_24.IMA  T2_47.IMA  T2_7.IMA   T2_92.IMA
T2_108.IMA T2_130.IMA  T2_153.IMA  T2_176.IMA  T2_25.IMA  T2_48.IMA  T2_70.IMA  T2_93.IMA
T2_109.IMA T2_131.IMA  T2_154.IMA  T2_177.IMA  T2_26.IMA  T2_49.IMA  T2_71.IMA  T2_94.IMA
T2_11.IMA  T2_132.IMA  T2_155.IMA  T2_178.IMA  T2_27.IMA  T2_5.IMA   T2_72.IMA  T2_95.IMA
T2_110.IMA T2_133.IMA  T2_156.IMA  T2_179.IMA  T2_28.IMA  T2_50.IMA  T2_73.IMA  T2_96.IMA
T2_111.IMA T2_134.IMA  T2_157.IMA  T2_18.IMA   T2_29.IMA  T2_51.IMA  T2_74.IMA  T2_97.IMA
T2_112.IMA T2_135.IMA  T2_158.IMA  T2_180.IMA  T2_3.IMA   T2_52.IMA  T2_75.IMA  T2_98.IMA
T2_113.IMA T2_136.IMA  T2_159.IMA  T2_181.IMA  T2_30.IMA  T2_53.IMA  T2_76.IMA  T2_99.IMA
T2_114.IMA T2_137.IMA  T2_16.IMA   T2_182.IMA  T2_31.IMA  T2_54.IMA  T2_77.IMA
T2_115.IMA T2_138.IMA  T2_160.IMA  T2_183.IMA  T2_32.IMA  T2_55.IMA  T2_78.IMA
T2_116.IMA T2_139.IMA  T2_161.IMA  T2_184.IMA  T2_33.IMA  T2_56.IMA  T2_79.IMA
T2_117.IMA T2_14.IMA   T2_162.IMA  T2_185.IMA  T2_34.IMA  T2_57.IMA  T2_8.IMA
T2_118.IMA T2_140.IMA  T2_163.IMA  T2_186.IMA  T2_35.IMA  T2_58.IMA  T2_80.IMA
T2_119.IMA T2_141.IMA  T2_164.IMA  T2_187.IMA  T2_36.IMA  T2_59.IMA  T2_81.IMA
```

We see that our 3-dimensional brain scan is composed many individual IMA files.

Reading Medical Images

Many of the functions we've used in the image processing toolbox have equivalents for DICOM and IMA files.

```
brain76Meta = dicominfo("T2_76.IMA") % Display meta data for DICOM/IMA image
```

```
brain76Meta = struct with fields:
```

```
    Filename: 'D:\Mina\Document\Coursera\Image Processing for Engineering and Science\Data\
    FileModDate: '04-Jan-2022 17:37:10'
    FileSize: 133222
    Format: 'DICOM'
    FormatVersion: 3
    Width: 256
    Height: 256
    BitDepth: 16
    ColorType: 'grayscale'
    FileMetaInformationGroupLength: 210
    FileMetaInformationVersion: [2x1 uint8]
    MediaStorageSOPClassUID: '1.2.840.10008.5.1.4.1.1.4'
    MediaStorageSOPInstanceUID: '1.3.6.1.4.1.9590.100.1.2.383125523312033013016153209983059058139'
    TransferSyntaxUID: '1.2.840.10008.1.2.1'
    ImplementationClassUID: '1.3.6.1.4.1.9590.100.1.3.100.9.4'
    ImplementationVersionName: 'MATLAB IPT 9.4'
    SpecificCharacterSet: 'ISO_IR 100'
    ImageType: 'ORIGINAL\PRIMARY\M\ND\NORM'
    InstanceCreationDate: ''
    InstanceCreationTime: ''
    SOPClassUID: '1.2.840.10008.5.1.4.1.1.4'
    SOPInstanceUID: '1.3.6.1.4.1.9590.100.1.2.383125523312033013016153209983059058139'
    StudyDate: '19390825'
    SeriesDate: '19390825'
    AcquisitionDate: '19390825'
    ContentDate: '19390825'
    StudyTime: '163228.526000'
    SeriesTime: '172035.762000'
```

AcquisitionTime: '165659.237500'
 ContentTime: '172035.920000'
 AccessionNumber: ''
 Modality: 'MR'
 Manufacturer: ''
 InstitutionName: 'Shiz University'
 InstitutionAddress: '1 Yellow Brick Road, Emerald City, Oz'
 ReferringPhysicianName: [1x1 struct]
 StationName: 'PoppyField'
 StudyDescription: 'IfOnlyaBrain'
 SeriesDescription: 'T2_tse_vfl_wave_1mm_R9_blip_RR'
 InstitutionalDepartmentName: 'Department'
 PerformingPhysicianName: [1x1 struct]
 OperatorsName: [1x1 struct]
 ManufacturerModelName: ''
 PatientName: [1x1 struct]
 PatientID: '99.99.99-99:99:99-DST-9.9.99.9.9999.9.9.99.99999'
 PatientBirthDate: '19040115'
 PatientSex: 'M'
 PatientAge: '107Y'
 PatientSize: 1.7907
 PatientWeight: 61.2350
 BodyPartExamined: 'BRAIN'
 ScanningSequence: 'SE'
 SequenceVariant: 'SK\SP'
 ScanOptions: 'PFP'
 MRAcquisitionType: '3D'
 SequenceName: 'spcR_260ns'
 AngioFlag: 'N'
 SliceThickness: 1
 RepetitionTime: 3200
 EchoTime: 406
 NumberOfAverages: 1
 ImagingFrequency: 123.2540
 ImagedNucleus: '1H'
 EchoNumbers: 1
 MagneticFieldStrength: 3
 NumberOfPhaseEncodingSteps: 223
 EchoTrainLength: 229
 PercentSampling: 100
 PercentPhaseFieldOfView: 100
 PixelBandwidth: 590
 DeviceSerialNumber: '45407'
 SoftwareVersions: ''
 ProtocolName: 'T2_tse_vfl_wave_1mm_R9_blip'
 TransmitCoilName: 'Body'
 AcquisitionMatrix: [4x1 uint16]
 InPlanePhaseEncodingDirection: 'ROW'
 FlipAngle: 120
 VariableFlipAngleFlag: 'Y'
 SAR: 0
 dBdt: 0
 PatientPosition: 'HFS'
 StudyInstanceUID: '1.3.12.2.1107.5.2.19.45407.30000018102520545318300000028'
 SeriesInstanceUID: '1.3.12.2.1107.5.2.19.45407.2018102617182732149201999.0.0.0'
 StudyID: '1'
 SeriesNumber: 21
 AcquisitionNumber: 1
 InstanceNumber: 76
 ImagePositionPatient: [3x1 double]
 ImageOrientationPatient: [6x1 double]
 FrameOfReferenceUID: '1.3.12.2.1107.5.2.19.45407.1.20181026163228673.0.0.4934'
 PositionReferenceIndicator: ''
 SliceLocation: -18.7376

```

        ImageComments: 'T2_tse_vfl_wave_1mm_R9_blip_RR+LW_181026_2/19900101/18.10.26-16:32:12-DST-
        SamplesPerPixel: 1
        PhotometricInterpretation: 'MONOCHROME2'
            Rows: 256
            Columns: 256
        PixelSpacing: [2x1 double]
        BitsAllocated: 16
        BitsStored: 16
        HighBit: 15
        PixelRepresentation: 0
        SmallestImagePixelValue: 0
        LargestImagePixelValue: 648
        WindowCenter: 210
        WindowWidth: 496
        PerformedProcedureStepStartDate: '19390825'
        PerformedProcedureStepStartTime: '163228.586000'
        PerformedProcedureStepID: 'MR19390825163228'
        PerformedProcedureStepDescription: 'FOLLOWYELLOWBRICKROAD'


```

```
imgBrain76 = dicomread("T2_76.IMA"); % Read in single image from 3D volume
```

Similar commands also exist for [Nifti](#) and [Analyze 7.5](#) files.

Once we've read in the medical image with the appropriate command, we can use the tools we've already used in this course. Use `imshow` to display `imgBrain76`:

```
imshow(uint8(imgBrain76))
```




If you're having trouble seeing the image, you may want to review [Viewing Meta-Data and non-uint8 Images](#).

3D Images as a Series of Images

Volumes are a series of images organized over x, y, and z spatial planes. One of the most intuitive ways to visualize this is to look at a series of 2D images taken at regular points over the third spatial plane.

Read in `T2_71.IMA` through `T2_75.IMA` and display them to get an idea of what the volume looks like in a section of the Y-Z (left to right, or "sagittal") plane.


```
imgBrain71 = dicomread("T2_71.IMA");  
imgBrain72 = dicomread("T2_72.IMA");  
imgBrain73 = dicomread("T2_73.IMA");  
imgBrain74 = dicomread("T2_74.IMA");  
imgBrain75 = dicomread("T2_75.IMA");  
  
montage({uint8(imgBrain71), uint8(imgBrain72), uint8(imgBrain73),  
uint8(imgBrain74), uint8(imgBrain75)})
```



Viewing volumes by looking at a series of images in one plane, or slices, is very popular with radiologists.

3D Images as Volumes

While viewing volumes as series of 2D images makes it easy to visualize 3D objects on the 2D surface of your screen, it also has drawbacks. Consider this image of a single slice of a stained tumor:



We can see the size of the removed tumor at a single cross-section, but usually we would want to see the whole tumor before it's removed in order to make treatment decisions.

Creating a Volume Object

Creating a volume allows us to see images on our screen with an approximation of the depth we see in the real world. It also allows us to calculate useful parameters such as volume and surface area.

As with our '/T2 MRI Scan' directory, 2D slices of a volume will usually be saved with sequential names. This enables MATLAB to assemble the images from a folder in the correct order to create a volume.

Before you can create the volume, you first need to set the path to the 'T2 MRI Scan' folder provided with the course files:

```
mrDir = uigetdir
```

```
mrDir =  
'D:\Mina\Document\Coursera\Image Processing for Engineering and Science\Data\MathWorks Images\T2 MRI Scan'
```

Now you are ready to load the volume:

```
vol = dicomreadVolume(fullfile(mrDir));
```

Removing extra dimensions

You might notice something strange about our new, 3D volume:

```
size(vol)
```

```
ans = 1x4  
    256    256     1    192
```

Why does our 3D volume have 4 dimensions? And why does our 3rd dimension have a length of 1? Because our images, like many medical images, are grayscale, the "color" dimension of our volume is empty, leading to a "singleton" dimension. Let's remove it and check the size again:

```
vol = squeeze(vol);  
size(vol)
```


```
ans = 1x3  
      256   256   192
```

That makes a lot more sense!

Viewing a Volume Object

Now we're finally ready to see our volume:

```
volshow(vol);
```



Now that is an impressive head. Let's check out some individual slices:

```
figure  
sliceViewer(vol)
```

```
ans =  
  sliceViewer with properties:  
  
    SliceDirection: [0 0 1]  
    SliceNumber: 97  
        Parent: [1x1 Panel]  
    Colormap: [256x3 double]  
  DisplayRange: [0 740]  
    ScaleFactors: [1 1 1]  
DisplayRangeInteraction: 'on'
```

You can scroll through slices from different planes using the `SliceDirection` parameter:

```
sliceViewer(vol,"SliceDirection","X")
```

```
ans =  
  sliceViewer with properties:  
    SliceDirection: 'X'  
    SliceNumber: 97  
    Parent: [1x1 Panel]  
    Colormap: [256x3 double]  
    DisplayRange: [0 740]  
    ScaleFactors: [1 1 1]  
    DisplayRangeInteraction: 'on'
```

Keep in mind, if your imported image isn't oriented correctly, the slice view direction might not correspond to the plane you expect. Try viewing the volume in different slice planes:

Once you have a volume object, you can view it, perform manual calculations, or import it from the workspace into the Volume Segmenter app.

Copyright 2022 The MathWorks, Inc.